

V2 INJECTOR

Skin Rejuvenation

Skin Hydration

Skin Elasticity



Discover the game changing Vital Injector2, a device that evenly delivers solutions to the skin using multi-needles. Various types of solutions can be injected to improve skin moisturization, elasticity, and texture. Engineered for excellence, it reduces pain and redefines precision, hygiene and operational efficiency.

The advanced Vital Injector2 boasts the capability to accommodate interchangeable syringes ranging from 1cc to 5cc. Its distinctive internal mechanism enables precise and consistent delivery of formulations, whether low or high in viscosity, in meticulously accurate dosages.

Specifically engineered to aid in administering stabilized hyaluronic acid filler, the Vital Injector2 plays a pivotal role in re-establishing skin hydro-balance. By enhancing skin elasticity and tone, it rejuvenates the skin, resulting in a revitalized and refreshed appearance. ARTG 433031.

Benefits:

- Increase skin moisture level
- Increase skin elasticity
- Increase skin texture and tone.
- Increase collagen production
- Rejuvenates the skin

- Revitalize and Refresh appearance
- Decrease pore size
- Decrease scar tissue
- Decrease sun damage

Three sessions are required initially, each a month apart, followed by a maintenance session recommended two to three times a year. Post treatment, your results should remain evident for at least 12 months before the effects of aging are noticeable again.

Results can be seen as soon as the skin has settled after the first treatment session. The skin should look noticeably healthier, less creased, better hydrated and more supple. With repeated sessions, your skin will look increasingly more radiant and feel firmer.

